QUICKIE®

Quickie® Chameleon / Chameleon HD
User Instruction Manual & Warranty

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I. FOREWORD

Dear Customer,

We are very happy that you have decided in favor of a high-quality product from SUNRISE MEDICAL.

This user's manual will provide numerous tips and ideas so that your new wheelchair can become a trustworthy and reliable partner in your life.

Maintaining close links with our customers is of great importance to us at Sunrise Medical. We therefore like to keep you up-to-date with our new and current developments. Keeping close to our customers also means fast service when you need replacement parts or accessories, or when you just have a question about your wheelchair– and with as little red tape as possible.

We want you to be satisfied with our products and service. Sunrise Medical therefore constantly works at continuous development of its products. For this reason, changes can occur in our range of products with regard to shape, technology, and fittings. Consequently, no claims can be construed from the data or pictures contained in this user's manual.

SUNRISE MEDICAL has been awarded the ISO 9001 Certificate, which affirms the quality of our products at every stage, from R & D to production.

Please contact your local, authorized SUNRISE MEDICAL dealer if you have any questions concerning the use, maintenance, or safety of your wheelchair.

In the case that there is no authorized dealer in your area or you have any questions, you can contact Sunrise Medical either in writing or by telephone at:

SUNRISE MEDICAL
Customer Service Department
7477 East Dry Creek Parkway
Longmont, CO 80503
800-333-4000

This manual gives information on features sold across different countries, without stating whether they are actually available in your country or are optional or standard features. For this information please refer to the order form/prescription form or your prescriber/supplier.

II. TABLE OF CONTENTS

I. FOREWORD

II. TABLE OF CONTENTS

III. WHEELCHAIR COMPONENTS

IV. NOTICE– READ BEFORE USE

V. GENERAL WARNINGS

VI. WARNINGS: FALLS & TIP-OVERS

VII. WARNINGS: FOR SAFE USE.....

VIII. WARNINGS: COMPONENTS & OPTIONS

IX. CHAIR BASICS.....

 A. Handling.....

 B. Quick-Release Axles for Rear Wheels

 C. Step Tubes

 D. Wheel Locks

 E. Footplates.....

 F. Legrests and Latch

 G. Casters

 H. Rear Axle

 I. Amputee Axle Position.....

 J. Backrest

 K. Armrest

 L. Seat Belt

 M. Anti-Tip Tubes.....

 N. Seat Depth.....

 O. Crutch Holder

 P. One-Arm Drive

 Q. Tires and Mounting.....

X. TROUBLESHOOTING.....

XI. MAINTENANCE

XII. TORQUE

XIII. WARRANTY

III. WHEELCHAIR COMPONENTS

We at SUNRISE MEDICAL have been awarded the ISO 9001 Certificate, which affirms the quality of our products at every stage, from R & D to production.

This product complies with the standards set forth in EU directives.

Optional equipment and accessories are available at extra charge.

1. Push handles
2. Backrest sling upholstery
3. Armrest
4. Seat sling upholstery—
shown w/optional cushion
5. Swing away Footrest
6. Casters
7. Caster housing
8. Fork
9. Quick-release axles – if so equipped
10. Wheel locks
11. Handrim
12. Rear wheel

QUICKIE CHAMELEON

SPECIFICATIONS:

	STANDARD	HEAVY DUTY
Weight	33 lbs (15 kg)	34 lbs (15.5kg)
Seat Widths	14, 15, 16, 17, 18, 19, 20	18, 19, 20, 22
Seat Depths	16, 17, 18	16, 17, 18, 19, 20
Front STFH– Std	16 - 20.5"	16 - 20.5"
Avail. Seat Angles	0, 3, 6 degrees	0, 3, 6 degrees
Frame Colors	Silver Texture, Black, Forest Green, Blue	Silver Texture, Black, Forest Green, Blue
Flip Back Armrests	Standard	Standard
Height Adj. Armrests	Optional	Optional
Rear Wheels	22 and 24"	22 and 24"
Front Casters	6 and 8"	6 and 8"
Forward Offset Axle	Standard	Standard
Amputee Axle	Optional	Optional
72° Swing-Away Hanger	Standard	Standard
Footrest Extension Tube Range	14.5" - 20"	14.5" - 20"
Elevating Legrests	Optional	Optional

IV. NOTICE– READ BEFORE USE

A. CHOOSE THE RIGHT CHAIR & SAFETY OPTIONS

Sunrise provides a choice of many wheelchair styles to meet the needs of the wheelchair rider. However, final selection of the type of wheelchair, options and adjustments rest solely with you and your health care advisor. Choosing the best chair and set-up for your safety depends on such things as:

1. Your disability, strength, balance and coordination.
2. The types of hazards you must overcome in daily use (where you live and work, and other places you are likely to use your chair).
3. The need for options for your safety and comfort (such as anti-tip tubes, positioning belts, or special seating systems).

B. REVIEW THIS MANUAL OFTEN

Before using this chair you, and each person who may assist you, should read this entire manual and make sure to follow all instructions. Review the warnings often, until they are second nature to you.

C. WARNINGS

The word “WARNING” refers to a hazard or unsafe practice that may cause severe injury or death to you or to other persons. The “Warnings” are in four main sections, as follows:

1. V. GENERAL WARNINGS

Here you will find a safety checklist and a summary of risks you need to be aware of before you ride this chair.

2. VI. WARNINGS: FALLS & TIP-OVERS

Here you will learn how to avoid a fall or tip-over while you perform daily activities in your chair.

3. VII. WARNINGS: FOR SAFE USE– ATTENDANTS

Here you will learn about practices for the safe use of your chair

4. VIII. WARNINGS: COMPONENTS & OPTIONS

Here you will learn about the components of your chair and options you can select for safety. Consult your authorized supplier and your health care advisor to help you choose the best set-up and options for safe use.

NOTE– Where they apply, you will also find “Warnings” in other sections of this manual.

Heed all warnings in this section. If you fail to do so a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

V. GENERAL WARNINGS

A. WEIGHT LIMIT

WARNING

NEVER exceed the weight limit of 275 pounds (125 kilograms) or 350 pounds (160 kilograms) for HD, for combined weight of rider and items carried. If you do exceed the limit, damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

B. WEIGHT TRAINING

WARNING

NEVER use this chair for weight training if total weight (rider plus equipment) exceeds 275 pounds (125 kilograms) or 350 pounds (160 kilograms) for HD. If you do exceed the limit, damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

C. GETTING TO KNOW YOUR CHAIR

WARNING

Every wheelchair is different. Take the time to learn the feel of this chair before you begin riding. Start slowly, with easy, smooth strokes. If you are used to a different chair, you may use too much force and tip over. If you use too much force, damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

D. TO REDUCE THE RISK OF AN ACCIDENT

WARNING

1. BEFORE riding, you should be trained in the safe use of this chair by your health care advisor.
2. Practice bending, reaching and transfers until you know the limit of your ability. Have someone help you until you know what can cause a fall or tip-over and how to avoid doing so.
3. Be aware that you must develop your own methods for safe use best suited to your level of function and ability.
4. NEVER try a new maneuver on your own until you are sure you can do it safely.
5. Get to know the areas where you plan to use your chair. Look for hazards and learn how to avoid them.
6. Use anti-tip tubes unless you are a skilled rider of this chair and are sure you are not at risk to tip over.
7. Propel your chair using the handrim only. Do not propel the chair by placing your hands on the rear tires.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

E. SAFETY CHECKLIST

WARNING

Before Each Use Of Your Chair:

1. Make sure the chair rolls easily and that all parts work smoothly. Check for noise, vibration, or a change in ease of use. (They may indicate low tire pressure, loose fasteners, or damage to your chair).
2. Repair any problem. Your authorized supplier can help you find and correct the problem.
3. Check to see that both quick-release rear axles are locked. When locked, the axle button will "pop out" fully. If not locked, the wheel may come off and cause you to fall.
4. If your chair has anti-tip tubes, lock them in place.

NOTE— See "Anti-Tip Tubes" (Section VIII) for times you should not use them.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

F. CHANGES & ADJUSTMENTS

WARNING

1. If you modify or adjust this chair, it may increase the risk of a tip-over UNLESS you make other changes as well.
2. Consult your authorized supplier BEFORE you modify or adjust your chair.
3. We recommend that you use anti-tip tubes until you adapt to the change, and are sure you are not at risk to tip over.
4. Unauthorized modifications or use of parts not supplied or approved by Sunrise may change the chair structure. This will void the warranty and may cause a safety hazard.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

G. ENVIRONMENTAL CONDITIONS

WARNING

1. Use extra care if you must ride your chair on a wet or slick surface. If you are in doubt, ask for help.
2. Contact with water or excess moisture may cause your chair to rust or corrode. This could cause your chair to fail.
 - a. Do not use your chair in a shower, pool or other body of water. The chair tubing and parts are not water-tight and may rust or corrode from the inside.
 - b. Avoid excess moisture (for example, do not leave your chair in a damp bathroom while taking a shower).
 - c. Dry your chair as soon as you can if it gets wet, or if you use water to clean it.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

H. TERRAIN

WARNING

1. Your chair is designed for use on firm, even surfaces such as concrete, asphalt and indoor floors and carpeting.
2. Do not operate your chair in sand, loose soil or over rough terrain. This may damage wheels or axles, or loosen fasteners of your chair.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

I. STREET USE

WARNING

In most states, wheelchairs are not legal for use on public roads. Be alert to the danger of motor vehicles on roads or in parking lots.

1. At night, or when lighting is poor, use reflective tape on your chair and clothing.
2. Due to your low position, it may be hard for drivers to see you. Make eye contact with drivers before you go forward. When in doubt, yield until you are sure it is safe.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

J. MOTOR VEHICLE SAFETY** WARNING**

Your wheelchair does not meet federal standards for motor vehicle seating.

1. NEVER let anyone sit in this chair while in a moving vehicle.
 - a. ALWAYS move the rider to an approved vehicle seat.
 - b. ALWAYS secure the rider with proper motor vehicle restraints.
2. In an accident or sudden stop the rider may be thrown from the chair. Wheelchair seat belts will not prevent this, and further injury may result from the belts or straps.
3. NEVER transport this chair in the front seat of a vehicle. It may shift and interfere with the driver.
4. ALWAYS secure this chair so that it cannot roll or shift.
5. Do not use any chair that has been involved in a motor vehicle accident.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

K. WHEN YOU NEED HELP** WARNING****For The Rider:**

Make sure that each person who helps you reads and follows all warnings and instructions that apply.

For Attendants:

1. Work with the rider's doctor, nurse or therapist to learn safe methods best suited to your abilities and those of the rider.
2. Tell the rider what you plan to do, and explain what you expect the rider to do. This will put the rider at ease and reduce the risk of an accident.
3. Make sure the chair has push handles. They provide secure points for you to hold the rear of the chair to prevent a fall or tip-over. Check to make sure push handle grips will not rotate or slip off.
4. To prevent injury to your back, use good posture and proper body mechanics. When you lift or support the rider or tilt the chair, bend your knees slightly and keep your back as upright and straight as you can.
5. Remind the rider to lean back when you tilt the chair backward.
6. When you descend a curb or single step, slowly lower the chair in one easy movement. Do not let the chair drop the last few inches to the ground. This may damage the chair or injure the rider.
7. To avoid tripping, unlock and rotate anti-tip tubes up, out of the way. ALWAYS lock the rear wheels and lock anti-tip tubes in place if you must leave the rider alone, even for a moment. This will reduce the risk of a tip over or loss of control of the chair.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

VI.WARNINGS: FALLS & TIP-OVERS**A. CENTER OF BALANCE****⚠ WARNING**

The point where this chair will tip forward, back, or to the side depends on its center of balance and stability. How your chair is set up, the options you select and the changes you make may affect the risk of a fall or tip-over.

1. The Most Important Adjustment Is:
The position of the rear wheels. The more you move the rear wheels forward, the more likely your chair will tip over backward.
2. The Center Of Balance Is Also Affected By:
 - a. A change in the set-up of your chair, including:
 - The distance between the rear wheels.
 - The amount of rear wheel camber.
 - The seat height and seat angle.
 - Backrest angle.
 - b. A change in your body position, posture or weight distribution.
 - c. Riding your chair on a ramp or slope.
 - d. The use of a back pack or other options, and the amount of added weight.
3. To Reduce The Risk Of An Accident:
 - a. Consult your doctor, nurse or therapist to find out what axle and caster position is best for you.
 - b. Consult your authorized supplier BEFORE you modify or adjust this chair. Be aware that you may need to make other changes to correct the center of balance.
 - c. Have someone help you until you know the balance points of your chair and how to avoid a tip-over.
 - d. Use anti-tip tubes.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

B. DRESSING OR CHANGING CLOTHES**⚠ WARNING**

Your weight may shift if you dress or change clothes while seated in this chair.

To reduce the risk of a fall or tip-over:

1. Rotate the front casters until they are as far forward as possible. This makes the chair more stable.
2. Lock anti-tip tubes in place. (If your chair does not have anti-tip tubes, back it up against a wall and lock both rear wheels).

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

C. WHEELIES**⚠ WARNING**

Doing a "wheelie" means balancing on the rear wheels of your chair while the front casters are in the air. It is dangerous to do a "wheelie" as a fall or tip-over may occur. However, if you do it safely, a "wheelie" can help you overcome curbs and obstacles.

1. Consult your doctor, nurse or therapist to find out if you are a good candidate to learn to do a "wheelie".
2. Do not attempt a "wheelie" UNLESS you are a skilled rider of this chair, or you have help.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

D. OBSTACLES**⚠ WARNING**

Obstacles and road hazards (such as potholes and broken pavement) can damage your chair and may cause a fall, tip-over or loss of control.

To avoid these risks:

1. Keep a lookout for danger - scan the area well ahead of your chair as you ride.
2. Make sure the floor areas where you live and work are level and free of obstacles.
3. Remove or cover threshold strips between rooms.
4. Install a ramp at entry or exit doors. Make sure there is not a drop off at the bottom of the ramp.
5. To Help Correct Your Center Of Balance:
 - a. Lean your upper body FORWARD slightly as you go UP over an obstacle.
 - b. Press your upper body BACKWARD as you go DOWN from a higher to a lower level.
6. If your chair has anti-tip tubes, lock them in place before you go UP over an obstacle.
7. Keep both of your hands on the handrims as you go over an obstacle.
8. Never push or pull on an object (such as furniture or a doorjamb) to propel your chair.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

E. REACHING OR LEANING**⚠ WARNING**

If you reach or lean it will affect the center of balance of your chair. This may cause you to fall or tip over. When in doubt, ask for help or use a device to extend your reach.

1. NEVER reach or lean if you must shift your weight sideways or rise up off the seat.
2. NEVER reach or lean if you must move forward in your seat to do so. Always keep your buttocks in contact with the backrest.
3. NEVER reach with both hands (you may not be able to catch yourself to prevent a fall if the chair tips).
4. NEVER reach or lean to the rear unless your chair has anti-tip tubes locked in place.
5. DO NOT reach or lean over the top of the seat back. This may damage one or both backrest tubes and cause you to fall.
6. If You Must Reach Or Lean:
 - a. Do not lock the rear wheels. This creates a tip point and makes a fall or tip-over more likely.
 - b. Do not put pressure on the footrests.

NOTE— Leaning forward puts pressure on the footrests and may cause the chair to tip if you lean too far.

- c. Move your chair as close as you can to the object you wish to reach.
- d. Do not try to pick up an object from the floor by reaching down between your knees. You are less likely to tip if you reach to the side of your chair.
- e. Rotate the front casters until they are as far forward as possible. This makes the chair more stable. Note: To do this: Move your chair past the object you want to reach, then back up alongside it. Backing up will rotate the casters forward.
- f. Firmly grasp a rear wheel or an armrest with one hand. This will help to prevent a fall if the chair tips.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

F. MOVING BACKWARD**⚠ WARNING**

Use extra care when you move your chair backward. Your chair is most stable when you propel yourself forward. You may lose control or tip over if one of the rear wheels hits an object and stops rolling.

1. Propel your chair slowly and smoothly.
2. If your chair has anti-tip tubes, make sure to lock them in place.
3. Stop often and check to be sure your path is clear.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

G. RAMPS, SLOPES & SIDEHILLS**⚠ WARNING**

Riding on a slope, which includes a ramp or sidehill, will change the center of balance of your chair. Your chair is less stable when it is at an angle. Anti-tip tubes may not prevent a fall or tip-over.

1. Do not use your chair on a slope steeper than 10%. (A 10% slope means: one foot in elevation for every ten feet of slope length.)
2. Always go as straight up and as straight down as you can. (Do not "cut the corner" on a slope or ramp.)
3. Do not turn or change direction on a slope.
4. Always stay in the CENTER of the ramp. Make sure ramp is wide enough that you are not at risk that a wheel may fall over the edge.
5. Do not stop on a steep slope. If you stop, you may lose control of your chair.
6. NEVER use rear wheel locks to try to slow or stop your chair. This is likely to cause your chair to veer out of control.
7. Beware Of:
 - a. Wet or slippery surfaces.
 - b. A change in grade on a slope (or a lip, bump or depression). These may cause a fall or tip-over.
 - c. A drop-off at the bottom of a slope. A drop-off of as small as 3/4 inch can stop a front caster and cause the chair to tip forward.
8. To Reduce The Risk Of A Fall Or Tip-Over:
 - a. Lean or press your body UPHILL. This will help adjust for the change in the center of balance caused by the slope or side hill.
 - b. Keep pressure on the handrims to control your speed on a down slope. If you go too fast you may lose control.
 - c. Ask for help any time you are in doubt.
9. Ramps At Home & Work - For your safety, ramps at home and work must meet all legal requirements for your area. We recommend:
 - a. Width. At least four feet wide.
 - b. guardrails. To reduce the risk of a fall, sides of ramp must have guardrails (or raised borders at least three inches high).
 - c. Slope. Not more than a 10% grade.
 - d. Surface. Flat and even, with a thin carpet or other non-skid material. (Make sure there is no lip, bump or depression.)
 - e. Bracing. Ramp must be STURDY. You may need bracing so ramp does not "bow" when you ride on it.
 - f. Avoid A Drop Off. You may need a section at the top or bottom to smooth out the transition.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

H. TRANSFER**⚠ WARNING**

It is dangerous to transfer on your own. It requires good balance and agility. Be aware that there is a point during every transfer when the wheelchair seat is not below you. To avoid a fall:

1. Work with your health care advisor to learn safe methods.
 - a. Learn how to position your body and how to support yourself during a transfer.
 - b. Have someone help you until you know how to do a safe transfer on your own.
2. Lock the rear wheels before you transfer. This keeps the rear wheels from rolling.

NOTE– This will NOT keep your chair from sliding away from you or tipping.

3. Make sure to keep pneumatic tires properly inflated. Low tire pressure may allow the rear wheel locks to slip.
4. Move your chair as close as you can to the seat you are transferring to. If possible, use a transfer board.
5. Rotate the front casters until they are as far forward as possible.
6. If you can, remove or swing footrests out of the way.
 - a. Make sure your feet do not catch in the space between the footrests.
 - b. Avoid putting weight on the footrests as this may cause the chair to tip.
7. Make sure armrests are out of the way and do not interfere.
8. Transfer as far back onto the seat surface as you can. This will reduce the risk that the chair will tip or move away from you.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

I. CURBS & STEPS**⚠ WARNING**

1. Each person who helps you should read and follow the warnings “For Safe Use” (Section VII).
2. Do not try to climb or descend a curb or step alone UNLESS you are a skilled rider of this chair and:
 - a. You can safely do a “wheelie”, and
 - b. You are sure you have the strength and balance to do so.
3. Unlock and rotate anti-tip tubes up, out of the way, so they do not interfere.
4. Do not try to climb a high curb or step (more than 4 inches high) UNLESS you have help. Doing so may cause your chair to exceed its balance point and tip over.
5. Go straight up and straight down a curb or step. If you climb or descend at an angle, a fall or tip-over is likely.
6. Return anti-tip tubes to locked position after overcoming obstacle.
7. Be aware that the impact of dropping down from a curb or step can damage your chair or loosen fasteners.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

J. STAIRS**⚠ WARNING**

1. NEVER use this chair on stairs UNLESS you have someone to help you. Doing so is likely to cause a fall or tip-over.
2. Persons who help you should read and follow the warnings “For Attendants” (Section V-K), and “For Safe Use” (Section VII).

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

K. ESCALATORS**⚠ WARNING**

NEVER use this chair on an escalator, even with an attendant. If you do, a fall or tip-over is likely.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

VII. WARNINGS: FOR SAFE USE– ATTENDANTS**⚠ WARNING**

Before you assist a rider, be sure to read the following warnings as well as the warnings “For Attendants” (section V-K) and follow all instructions that apply. Be aware that you will need to learn safe methods best suited to your abilities.

A. LEARNING TO DO A “WHEELIE”**⚠ WARNING**

Follow these steps to help the rider learn to do a “wheelie”:

1. Read and follow the warnings for “WHEELIES”. (See section VI-C)
2. Stand at the rear of the chair. You must be able to move with the chair to prevent a tip-over.
3. Keep your hands BENEATH the push handles, ready to catch the rider if the chair tips back.
4. When you are ready, have the rider place his or her hands on the forward part of the handrims.
5. Have the rider make a quick BACKWARD movement of the rear wheels, quickly followed by a hard FORWARD thrust. (This will tilt the chair up on the balance point of the rear wheels.)
6. Have the rider make small movements of the handrims to stay in balance.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

B. DESCENDING A CURB OR SINGLE STEP**⚠ WARNING**

Follow these steps to help a rider descend a curb or single step going BACKWARD:

1. Stay at the rear of the chair.
2. Several feet before you reach the edge of the curb or step, turn the chair around and pull it backward.
3. While looking over your shoulder, carefully step back until you are off the curb or stair and standing on the lower level.
4. Pull the chair toward you until the rear wheels reach the edge of the curb or step. Then allow the rear wheels to slowly roll down onto the lower level.
5. When the rear wheels are safely on the lower level, tilt the chair back to its balance point. This will lift the front casters off the curb or step.
6. Keep the chair in balance and take small steps backward. Turn the chair around and gently lower front casters to the ground.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

C. CLIMBING A CURB OR SINGLE STEP**⚠ WARNING**

Follow these steps to help the rider climb a curb or single step going FORWARD:

1. Stay behind the chair.
2. Face the curb and tilt the chair up on the rear wheels so that the front casters clear the curb or step.
3. Move forward, placing the front casters on the upper level as soon as you are sure they are past the edge.
4. Continue forward until the rear wheels contact the face of the curb or step. Lift and roll the rear wheels to the upper level.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

D. CLIMBING STAIRS** WARNING**

1. Use at least two attendants to move a chair and rider up stairs.
2. Move the chair and rider BACKWARD up the stairs.
3. The person at the rear is in control. He or she tilts the chair back to its balance point.
4. A second attendant at the front firmly grasps a non-detachable part of the front frame and lifts the chair up and over one stair at a time.
5. The attendants move to the next stair up. Repeat for each stair, until you reach the landing.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

E. DESCENDING STAIRS** WARNING**

1. Use at least two attendants to move a chair and rider down stairs.
2. Move the chair and rider FORWARD down the stairs.
3. The person at the rear is in control. He or she tilts the chair to the balance point of the rear wheels and rolls it to the edge of the top step.
4. A second attendant stands on the third step from the top and grasps the chair frame. He or she lowers the chair one step at a time by letting the rear wheels roll over the stair edge.
5. The attendants move to the next stair down. Repeat for each stair, until you reach the landing.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

VIII. WARNINGS: COMPONENTS & OPTIONS**A. ANTI-TIP TUBES (OPTIONAL)****⚠ WARNING**

Anti-tip tubes can help keep your chair from tipping over backward in most normal conditions.

1. Sunrise Recommends Use Of Anti-Tip Tubes:
 - a. UNLESS you are a skilled rider of this chair and are sure you are not at risk to tip over.
 - b. Each time you modify or adjust your chair. The change may make it easier to tip backward. Use anti-tip tubes until you adapt to the change, and are sure you are not at risk to tip over.
2. When locked in place (in the “down” position) anti-tip tubes should be BETWEEN 1 to 2 inches (3cm to 5cm) off the ground.
 - a. If set too HIGH, they may not prevent a tip-over.
 - b. If set too LOW, they may “hang up” on obstacles you can expect in normal use. If this occurs, you may fall or your chair may tip over.
3. Keep Anti-Tip Tubes Locked In Place UNLESS:
 - a. You have an attendant, or
 - b. You have to climb or descend a curb, or overcome an obstacle, and can safely do so without them. At these times, make sure anti-tip tubes are up, out of the way.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

B. ARMRESTS**⚠ WARNING**

Armrests detach and will not bear the weight of this chair.

1. NEVER lift this chair by its armrests. They may come loose or break.
2. Lift this chair only by non-detachable parts of the main frame.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

C. CUSHIONS & SLING SEATS**⚠ WARNING**

1. Sling seats and standard foam cushions are not designed for the relief of pressure.
2. If you suffer from pressure sores, or if you are at risk that they will occur, you may need a special seat system or a device to control your posture. Consult your doctor, nurse or therapist to find out if you need such a device for your well-being.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

D. FASTENERS**⚠ WARNING**

Many of the screws, bolts and nuts on this chair are special high-strength fasteners. Use of improper fasteners may cause your chair to fail.

1. ONLY use fasteners provided by an authorized supplier (or ones of the same type and strength, as indicated by the markings on the heads).
2. Over- or under-tightened fasteners may fail or cause damage to chair parts.
3. If bolts or screws become loose, tighten them as soon as you can.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

E. FOOTRESTS**⚠ WARNING**

1. At the lowest point, footrests should be AT LEAST 2 1/2 INCHES (64 mm) off the ground. If set too LOW, they may “hang up” on obstacles you can expect to find in normal use. This may cause the chair to stop suddenly and tip forward.
2. To Avoid A Trip Or Fall When You Transfer:
 - a. Make sure your feet do not “hang up” or get caught in the space between the footrests.
 - b. Avoid putting weight on the footrests, as the chair may tip forward.
3. NEVER lift this chair by the footrests. Footrests detach and will not bear the weight of this chair. Lift this chair only by non-detachable parts of the main frame.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

F. PNEUMATIC TIRES**⚠ WARNING**

Proper inflation extends the life of your tires and makes your chair easier to use.

1. Do not use this chair if any of the tires is under- or over-inflated. Check weekly for proper inflation level, as listed on the tire sidewall.
2. Low pressure in a rear tire may cause the wheel lock on that side to slip and allow the wheel to turn when you do not expect it.
3. Low pressure in any of the tires may cause the chair to veer to one side and result in a loss of control.
4. Over-inflated tires may burst.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

G. POSITIONING BELTS (OPTIONAL)**⚠ WARNING**

Use positioning belts ONLY to help support the rider’s posture. Improper use of these belts may cause severe injury to or death of the rider.

1. Make sure the rider does not slide down in the wheelchair seat. If this occurs, the rider may suffer chest compression or suffocate due to pressure from the belts.
2. The belts must be snug, but must not be so tight that they interfere with breathing. You should be able to slide your open hand, flat, between the belt and the rider.
3. A pelvic wedge or a similar device can help keep the rider from sliding down in the seat. Consult with the rider’s doctor, nurse or therapist to find out if the rider needs such a device.
4. Use positioning belts only with a rider who can cooperate. Make sure the rider can easily remove the belts in an emergency.
5. NEVER Use Positioning Belts:
 - a. As a patient restraint. A restraint requires a doctor’s order.
 - b. On a rider who is comatose or agitated.
 - c. As a motor vehicle restraint. In an accident or sudden stop the rider may be thrown from the chair. Wheelchair seat belts will not prevent this, and further injury may result from the belts or straps.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

H. POWER DRIVE**⚠ WARNING**

Do not install a power drive on your wheelchair. If you do:

1. It will affect the center of balance of your chair and may cause a fall or tip-over.
2. It will alter the frame and void the warranty.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

I. PUSH HANDLES (OPTIONAL)**⚠ WARNING**

When you have an attendant, make sure that this chair has push handles.

1. Push handles provide secure points for an attendant to hold the rear of this chair, to prevent a fall or tip-over. Make sure to use push handles when you have an attendant.
2. Check to make sure push handle grips will not rotate or slip off.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

J. QUICK-RELEASE AXLES (OPTIONAL)**⚠ WARNING**

1. Do not use this chair UNLESS you are sure that both quick-release rear axles are locked. An unlocked axle may come off during use and cause a fall.
2. An axle is not locked until the quick-release button pops out fully.
An unlocked axle may come off during use, resulting in a fall, tip-over or loss of control and cause severe injury to the rider or others.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

K. REAR WHEELS**⚠ WARNING**

A change in set-up of the rear wheels will affect the center of balance of your chair.

1. The farther you move the rear axles FORWARD, the more likely it is that your chair will tip over backward.
2. Consult your doctor, nurse or therapist to find the best rear axle set-up for your chair.
Do not change the set-up UNLESS you are sure you are not at risk to tip over.
3. Adjust the rear wheel locks after you make any change to the rear axles.
 - a. If you fail to do so, the locks may not work.
 - b. Make sure lock arms embed in tires at least 1/8 inch (3 mm) when locked.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

L. REAR WHEEL LOCKS**⚠ WARNING**

Rear wheel locks are NOT designed to slow or stop a moving wheelchair. Use the only to keep the rear wheels from rolling when your chair is at a complete stop.

1. NEVER use rear wheel locks to try to slow or stop your chair when it is moving. Doing so may cause you to veer out of control.
2. To keep the rear wheels from rolling, always set both rear wheel locks when you transfer to or from your chair.
3. Low pressure in a rear tire may cause the wheel lock on that side to slip and may allow the wheel to turn when you do not expect it.
4. Make sure lock arms embed in tires at least 1/8 inch (3 mm) when locked. If you fail to do so, the locks may not work.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

M. MODIFIED SEAT SYSTEMS** WARNING**

Use of a seat system not approved by Sunrise may alter the center of balance of this chair. This may cause the chair to tip over.

1. Do not change the seat system of your chair UNLESS you consult your authorized supplier first.
2. Use of a seat system not approved by Sunrise may affect the folding mechanism of this chair.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

N. UPHOLSTERY FABRIC** WARNING**

1. Replace worn or torn fabric of seat and seat back as soon as you can. If you fail to do so, the seat may fail.
2. Sling fabric will weaken with age and use. Look for fraying or thin spots, or stretching of fabrics at mounting holes.
3. "Dropping down" into your chair will weaken fabric and result in the need to inspect and replace the seat more often.
4. Be aware that laundering or excess moisture will reduce flame retardation of the fabric.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

IX. CHAIR BASICS

A. HANDLING

Folding up (Fig. 1)

First remove the seat cushion from the wheelchair and flip up the individual footplates.

Take hold of the sling from underneath in the middle and pull upwards. The wheelchair will then fold together.

To fold your wheelchair so that it is as small as possible, e.g., to put into a car, you can remove the footrests. For this purpose, press the latch and swing the footrest to the side. Then the footrest only has to be drawn out of the frame tubing.

Unfolding (Fig. 2)

Press down on the seat rails (see fig. 2). Your wheelchair will then unfold. Now snap the seat rails into position in the seat saddles. This can be facilitated by slightly tipping your wheelchair to the side. In this way, the weight can be taken off one rear wheel so that the wheelchair can be more easily unfolded.

Be careful not to get your fingers caught in the cross-tube assembly.

B. QUICK-RELEASE AXLES FOR REAR WHEELS – IF SO EQUIPPED (FIG. 3)

The rear wheels (24", 22") may be equipped with quick-release axles. The wheels can be installed or removed without using tools. To remove a wheel, simply depress the quick-release button on the axle (A) and pull the wheel out.

⚠ WARNING

Hold the quick-release button on the axle depressed when inserting the axle into the frame to mount the rear wheels. Release the button to lock the wheel in place. The quick-release button (A) must snap back to its original position to ensure proper engagement.

C. STEP TUBES (FIG. 4)

Step tubes are used by attendants to tip a wheelchair over an obstacle. Simply step on the tube to push a wheelchair, for example, over a curb or step.

D. WHEEL LOCKS (FIG. 5)

Your wheelchair is equipped with two wheel locks. They are applied directly against the tires.

To engage, press both wheel-lock levers (B) forward against the stops. To release the wheel, pull the levers back to their rearward positions. Note: For "pull-to-lock" wheel locks, pull back on both levers to engage and push forward to release.

Holding power of wheel locks will decrease with:



- Worn tire tread
- Pneumatic tire pressure that is too low
- Wet tires
- Improperly adjusted wheel locks.

The wheel locks have not been designed to be used as brakes for a moving wheelchair. The wheel locks should therefore never be used to brake a moving wheelchair. Always use the handrims for braking.

Make sure that the interval between the tires and wheel locks is adjusted to provide sufficient holding force. To readjust, loosen clamp screw found on underside of mechanism. Set the appropriate interval. Tighten screw. (See section XII. TORQUE) Note: a T30 Torx screwdriver is required for wheel lock adjustment.

⚠ WARNING

After each adjustment of the rear wheels, check the interval to the wheel locks and readjust if necessary.

Unilateral Wheel Lock– if so equipped

Your chair may be equipped with a unilateral, push-to-lock wheel lock mechanism. This mechanism allows engagement of both rear wheel locks by pushing forward on a single lever. Pulling rearward on the single lever releases both wheel locks.

Extension for the Wheel Lock Lever– if so equipped (Fig. 6)

The extension for the wheel lock lever can be removed or folded down. The longer lever helps to minimize the effort needed to set the wheel locks.

⚠ WARNING

Mounting the wheel lock too close toward the wheel will result in a higher effort to operate. This might cause the wheel lock lever extension to break!

⚠ WARNING

Leaning onto the wheel lock lever extension while transferring will cause the lever to break!

⚠ WARNING

Splashing water from tires might cause the wheel lock to malfunction.

E. FOOTPLATES (Fig. 7)

The footplates can be flipped up to facilitate getting in and getting out of your wheelchair.

Lower Leg Length Adjustment

By loosening the cap screws (A) on both footrest hangers, the legrests can be adjusted to fit any lower leg length. Loosen the screws, set the tubing with legrest in the desired position, and then tighten the screws. (See section XII. TORQUE).

⚠ WARNING

A minimum interval of 2 1/2 inches (64mm) from the ground should always be maintained to provide sufficient clearance. If set too low, they may “hang up” on obstacles. This may cause the chair to stop suddenly and tip forward.

Angle Adjustable Footplates – if so equipped

Your chair may be equipped with angle adjustable footplates. To adjust, loosen fasteners securing footplate sufficient to allow rotation. Rotate to appropriate angle and re-tighten fasteners.

F. LEGRESTS AND LATCH (Fig. 8)

The legrest can be swung in, underneath the seat sling, or outward to remove.




Mount legrest by fully inserting hanger plug (B) into frame with the footplate(s) oriented toward the in/outside. Then rotate the legrest toward the midline until it locks into place.

Test to make sure that the legrest is properly locked into place.

To remove the legrest, engage the lever (C) rearward to disengage latch, then swing the footrest toward the in/outside. Lift to remove it.

⚠ WARNING

The legrests should not be used to lift or to carry the wheelchair with an occupant. Lift chair only by non-detachable parts of the main frame.

Elevating Legrest– if so equipped (Fig. 9)**To elevate:**

Pull the footrest upward to adjust it to the desired height. Footrest will engage automatically in place.

To lower:

Lift the lower leg and press the release lever (A) to unlock the footrest. You may now lower the footrest. As soon as you let go of the lever, the footrest will lock back into place.

To Remove:

Depress release lever (B), swing out and then lift off.

⚠ WARNING

Keep hands clear of the adjustment mechanism between the frame and the movable parts of the footrest while elevating or lowering the footrest.

⚠ WARNING

Legrests are not to be used for lifting or carrying the wheelchair with an occupant.

Stump Pad– if so equipped

Your chair may be equipped with an amputee stump pad. Angle support appropriately for maximum comfort and support.

G. CASTERS (Figs. 10, 11, 12)**⚠ WARNING**

Casters should always be adjusted by an authorized dealer.

Caster Angle Adjustment

The angle of the caster should be checked whenever there is an adjustment to the seat to floor height (STFH) to ensure proper orientation.

To adjust, remove the plastic cap and caster nut (C). Remove stem bolt from caster housing. Remove hex collar (D) and reorient as follows:

Level seat (i.e. front and rear STFH are equal) = Double notch on top and to rear (E)

Rear STFH is 1" lower than front STFH = Single notch on top and to rear (F)

Rear STFH is 2" lower than front STFH = Double notch on top and to front (G)

To ensure proper alignment desired notch on hex collar must be on the top of the housing and on the rearward/forward edge as indicated above on both sides of chair. Replace stem bolt and caster nut to appropriate torque setting (see section XII. TORQUE).

Caster Position – Front Seat to Floor Height Adjustment – see Seat Height Matrix.

Front STFH is adjusted by:




1. Using a smaller/larger diameter caster wheel (6" or 8").
2. Remove caster axle (H) and position in alternate fork position (1").
3. Reorient caster stem spacers (I).

Remove plastic cap and caster nut (C). Remove stem bolt from caster housing (J).

Reorient 1" spacers for desired STFH. Options are:

- 2 spacers above housing = lowest
- 1 spacer above and 1 below housing = middle
- 2 spacers below housing = highest

NOTE– Alteration to front STFH, unless done in a 1 to 1 ratio with the rear STFH, will require Caster Angle Adjustment as previously described.

SEAT HEIGHT MATRIX (FIGS. 13, 14, 15)

Front STFH	Rear STFH	Rear Wheel	Rear Axle Position (Fig. 14)	Caster Wheel	Caster Collar, Angle Setting (Fig. 13)	Caster Axle Position (Fig. 15)	Caster Spacers Below Housing (Fig. 15)
16.5	14.5	20"	1	6"	G	1	0
16.5	15.5	22"	1	6"	F	1	0
16.5	16.5	22"	2	6"	E	1	0
17.5	15.5	22"	1	6"	G	2	0
17.5	16.5	24"	1	6"	F	2	0
17.5	17.5	24"	2	6"	E	2	0
18.5	16.5	24"	1	6"	G	1	2
18.5	17.5	24"	2	6"	F	1	2
18.5	18.5	24"	3	6"	E	1	2
19.5	17.5	24"	2	6"	G	2	2
19.5	18.5	24"	3	6"	F	2	2
19.5	19.5	24"	4	6"	E	2	2
18.5	16.5	24"	1	8"	G	2	0
18.5	17.5	24"	2	8"	F	2	0
18.5	18.5	24"	3	8"	E	2	0
19.5	17.5	24"	2	8"	G	2	1
19.5	18.5	24"	3	8"	F	2	1
19.5	19.5	24"	4	8"	E	2	1
20.5	18.5	24"	3	8"	G	2	2
20.5	19.5	24"	4	8"	F	2	2

12" Rear Wheels

16.5	16.5	12"	T	6"	E	1	0
17.5	16.5	12"	T	6"	F	2	0
18.5	16.5	12"	T	6"	G	1	2
18.5	16.5	12"	T	8"	G	2	0

Ultra Hemi Bracket

14.5	14.5	20"	1	6"	E	1	0
15.5	14.5	20"	1	6"	F	2	0
15.5	15.5	22"	1	6"	E	2	0

All measurements are approximate. Shaded areas are factory available configurations.

CAUTION— Measured to seat sling, without seat cushion.

H. REAR AXLE (Fig. 16)

Axle Position - Rear Seat to Floor Height Adjustment – see Seat Height Matrix.

Rear STFH is adjusted by:

- Using a larger/smaller diameter rear wheel (22" or 24").
- Repositioning the axle receiver in the axle plate.

Remove rear wheel axle. Remove axle receiver (A) from current position and replace in desired position using appropriate torque setting (see section XII. TORQUE).

NOTE— Alteration to rear STFH, unless done in a 1 to 1 ratio with the front STFH, will require Caster Angle Adjustment as previously described.

⚠ WARNING

Wheel locks must be adjusted whenever rear axle position is altered.






I. AMPUTEE AXLE POSITION (Fig. 17)

The longer the wheel base, the greater the stability of the wheelchair. Reposition the axle plate to the rear (B) in order to extend the wheel base.

When you have completed the adjustment, make sure that all screws have been properly tightened (See section XII. TORQUE).

⚠ WARNING

Wheel locks must be adjusted for the new position.

J. BACKREST

Height adjustable backrest (Fig. 18)

The backrest height can be adjusted to 5 different positions (38 – 48 cm, 15" –19"). Remove the bolt (A) and move the back tube into the desired position. Replace and tighten bolt to appropriate torque setting. (See section XII. TORQUE).

Angle-Adjustable Back Posts– if so equipped

The backrest can be adjusted to several different angles. A nut is located on the backtube, allowing you to set the angle (See section XII. TORQUE).

Angle adjustment of the back post will require adjustment of the armrest “hook” support.

For your own safety, we recommend that adjustments only be made by authorized dealers.

⚠ WARNING

Anti-tip tubes should be used in conjunction with angle-adjustable back posts.

Reclining Backrest– if so equipped


Your chair may be equipped with a reclining backrest. An attendant may adjust the angle of the backrest by engaging and releasing the levers.

⚠ WARNING

Reclining backrest should be used in conjunction with anti-tip tubes.

Tension Adjustable Backrest Upholstery– if so equipped

The sling can be adapted to a desired degree of tension by adjusting the Velcro fasteners.

K. ARMREST

Standard armrest, Flip-Up, Detachable with Desk or Full Armrests (Fig. 19)

To flip the armrest up, push the lever (B) to unlock and lift. Armrest will stand up behind backrest.

The height of the armpad can be set in an alternate position by reconfiguring the spacers. To do so, you must first loosen the screws accessed through bottom of arm body, reconfigure the spacers and then retighten the screws. (See section XII. TORQUE).

⚠ WARNING

Neither the armrests nor the armpads are to be used for lifting or carrying the wheelchair.

Height adjustable armrest, Flip-Up, Detachable with Desk or Full Armrests (Fig. 19)

To flip the armrest up, press the lever (B) to unlock and lift.

Adjust the height of the armrest as follows:

Press the button found on the underside of the armpad and move the armpad into the desired position. Release the button and push the armpad down until you can hear it click into place. Always make sure that the armrests have locked into place properly.

⚠ WARNING

Neither the armrest nor the armpads are to be used for lifting or carrying the wheelchair.

L. SEAT BELT

The seat belt is a hip belt that enhances wheelchair safety and proper positioning. The seat belt is fastened to the frame as shown in the picture.

NOTE– If an after-market belt is to be used, convenient attachment points may be purchased. Contact Sunrise Medical Customer Service and ask for “Positioning Belt Slide Clamps”.

M. ANTI-TIP TUBES

Anti-tip tubes provide additional safety for inexperienced users when they are still learning how to operate their wheelchairs.

Anti-tip tubes prevent a wheelchair from tipping over backward.

Pressing the release button allows the anti-tip tubes to be rotated upward or removed. An interval of 3 cm to 5 cm (1 to 2 inches) should be maintained between the anti-tip tubes and the ground.



You must rotate the anti-tip tubes upwards when going up and down large obstacles (such as a curb) to prevent them from touching the ground. Rotate anti-tip tubes down into functional position after overcoming obstacle.

N. SEAT DEPTH (Figs. 20, 21, 22)

Seat depth is adjusted in two ways:

1. By removing the clips (A) and loosening the cross-tube fasteners (B), the cross-tube assembly (C) can be moved along the frame, thereby allowing the seat depth to be changed (depending on the position of the back tubes). Make sure that the clips (A) always snap into the holes provided in the frame.
2. The seat depth also can be adjusted via the rear frame assembly in one inch increments. Remove the wheels and armrests and check if the cross brace is in the desired position (slide cross brace if necessary as described above). Remove the axle plate screws (D). Move the rear frame assembly into the desired position and replace and tighten all 3 screws. Finally, move the front receiver of the armrest into the required position.

A 2" long seat rail plug (E) and longer seat sling are available for 18" seat depth (18" & 20" on HD).

⚠ WARNING

Wheel locks must be adjusted for the new position.

O. CRUTCH HOLDER– OPTIONAL

This device permits crutches to be transported directly on a wheelchair. It has a Velcro® loop to fasten crutches or other aids.

P. ONE-ARM DRIVE– IF SO EQUIPPED

Your chair may be equipped with a One-Arm Drive unit. Rotation of both handrims on one side of the chair will facilitate independent propulsion.

Q. PNEUMATIC TIRES AND MOUNTING

When equipped with pneumatic tires, always make sure you that you maintain the correct tire pressure, as this can have an effect on wheelchair performance. If the tire pressure is too low, rolling resistance will increase, requiring more effort to move the chair forward; low tire pressure also has a negative impact on maneuverability. If the tire pressure is too high, the tire could burst.

The correct pressure for a given tire is printed on the surface of the tire itself.

Tires can be mounted the same way as an ordinary bicycle tire.

Before installing a new inner tube, you should always make sure that the base of the rim and the interior of the tire are free of foreign objects.

Check the pressure after mounting or repairing a tire. It is critical to your safety and to the wheelchair's performance that regulation air pressure be maintained and that tires be in good condition.

X. TROUBLESHOOTING

SYMPTOMS							SOLUTIONS
	Left Turn In Chair	Right Turn In Chair	Looseness In Chair	Sluggish Turning	Squeaks And Rattles	Caster Flutter	
✓	✓	✓	✓				Make sure tire pressure is correct and equal in both rear tires and front caster tires, if pneumatic.
		✓	✓	✓	✓		Make sure all nuts and bolts are snug.
		✓		✓			Make sure all spokes and nipples are tight on radial spoke wheels.
				✓			Use Tri-flow Lubricant (Teflon®-based) between all modular frame connections and parts.
✓	✓					✓	Make sure both front casters touch the ground simultaneously. If they do not, add the proper spacers between the bottom bearing of the caster plate and fork stem nut until they do. Make sure you check for this problem on a flat surface.

XI. MAINTENANCE

A. INTRODUCTION

1. Proper maintenance will improve performance and extend the useful life of your chair.
2. Clean your chair regularly. This will help you find loose or worn parts and make your chair easier to use.
3. Repair or replace loose, worn, bent or damaged parts before using this chair.
4. To protect your investment, have all major maintenance and repair work done by your authorized supplier.

WARNING

If you fail to heed these warnings your chair may fail and cause severe injury to the rider or others.

1. Inspect and maintain this chair strictly per chart.
2. If you detect a problem, make sure to service or repair the chair before use.
3. At least once a year, have a complete inspection, safety check and service of your chair made by an authorized supplier.

B. MAINTENANCE CHART

Weekly	3 Months	6 Months	Annually	CHECK...
✓				Tire inflation level
✓				Wheel locks
	✓			All fasteners for wear and tightness
		✓		Armrests
		✓		Axle and axle sleeves
	✓			Wheels, tires and spokes
	✓			Casters
	✓			Anti-tip tubes
		✓		Frame
	✓			Upholstery
			✓	Service by authorized supplier

C. MAINTENANCE TIPS

1. Axles & Axle Sleeves:

Check axles and axle sleeves every six months to make sure they are tight. Loose sleeves will damage the axle plate and will affect performance.

2. Tire Air Pressure:

Check air pressure in pneumatic tires at least ONCE A WEEK. The wheel locks will not grip properly if you fail to maintain the air pressure shown on tire sidewall.

D. CLEANING**1. Paint Finish**

- a. Clean the paint finish with mild soap at least once a month.
- b. Protect the paint with a coat of non-abrasive auto wax every three months.

2. Axles and Moving Parts

- a. Clean around axles and moving parts WEEKLY with a slightly damp (not wet) cloth.
- b. Wipe off or blow away any fluff, dust, or dirt on axles or moving parts.

NOTE– You do not need to grease or oil the chair.

3. Upholstery

- a. Hand-wash (machine washing may damage fabric).
- b. Drip-dry only. DO NOT machine dry as heat will damage fabric.

E. STORAGE TIPS

1. Store your chair in a clean, dry area. If you fail to do so, parts may rust or corrode.
2. Before using your chair, make sure it is in proper working order. Inspect and service all items on the “Maintenance Chart”.
3. If stored for more than three months, have your chair inspected by an authorized supplier before use.

F. NECESSARY TOOLS

- 19mm socket
- 10mm wrench or socket
- T30 Torx screwdriver/wrench
- 6mm hex key
- 5mm hex key
- 4mm hex key
- Large Adjustable Wrench
- Phillips screwdriver
- Torque wrench– recommended

XII. TORQUE

Sunrise recommends the use of a torque wrench when adjusting and maintaining this chair.

NOTE– A torque setting is the optimum tightening which should be made on a particular fastener. A torque wrench that measures in in-lbs, ft-lbs, or Nm should be used when tightening fasteners on your chair.

RECOMMENDED TORQUE SETTINGS:

	in-lbs	ft-lbs	Nm
Caster Bracket Bolts	88	7.4	10
Molded Parts	44	3.7	5
Caster Stem Bolt	531	44	60
Threaded Axle Sleeve	265	22	30
Normal Steel Parts	62	5.2	7

XIII. WARRANTY**A. LIFETIME LIMITED WARRANTY**

Sunrise warrants the frame and quick-release axles of this wheelchair against defects in materials and workmanship for the life of the original purchaser.

B. FOR TWO (2) YEARS

We warrant all Sunrise-made parts and components of this wheelchair against defects in materials and workmanship for two years from the date of first consumer purchase.

C. LIMITATIONS

1. We do not warrant:
 - Tires and tubes, upholstery, armrests and push-handle grips.
 - Damage from neglect, misuse, or from improper installation or repair.
 - Damage from exceeding weight limit of 275 pounds (125 kg) or 350 pounds (160 kg) for HD.
2. This warranty is VOID if the original chair serial number tag is removed or altered.
3. This warranty applies in the USA only. Check with your supplier to find out if international warranties apply.

D. WHAT WE WILL DO

Our sole liability is to repair or replace covered parts. This is the exclusive remedy for consequential damages.

E. WHAT YOU MUST DO

1. Obtain from us through an authorized supplier, while this warranty is in effect, prior approval for return or repair of covered parts.
2. Return the wheelchair or part(s), freight pre-paid, to Sunrise Medical at: 33554 Pin Oak Parkway, Avon Lake, OH 44012-2320.
3. Pay the cost of labor to repair, remove or install parts.

F. NOTICE TO CONSUMER

1. If allowed by law, this warranty is in place of any other warranty (written or oral, express or implied, including a warranty of merchantability, or fitness for a particular purpose).
2. This warranty gives you certain legal rights. You may also have other rights that vary from state to state.

Record your serial number here for future reference:

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Sunrise Medical, Inc.