



"Experience the performance and durability that I demand."

Bryan Anderson
Quantum Rehab Spokesperson



The Litestream[®] XF ultralight manual wheelchair uses a unique, innovative rigidizing double X folding frame that delivers the versatility and performance of a solid frame while folding to a width of only 9" for easy transport. The Litestream XF accepts a broad selection of accessories to meet a wide variety of needs and personal preferences while remaining an excellent overall value.



(US) 800-800-8586 • (Canada) 888-570-1113 Find out more at www.quantumrehab.com



SPECIFICATIONS	
SADMERC Code:	K0005
Product Weight:	27 lbs.
Weight Capacity:	250 lbs.
Seat to Floor Heights:	15"-19.5" front 14.5"-19" rear
Seat Sizes:	Width: 14"-21" Depth: 14"-20" (Available in 1" increments)







OPTIONS	
Rear Wheels:	24", 26" Spinergy wire spoke 20", 22", 24" Composite mag 24" Spinergy SPOX (available in red, blue, yellow, black & white) ¹
Tires:	100 PSI Primo tire Aggressive knobby tire Everyday gray tire Polyurethane tire
Handrims:	Aluminum Vinyl-coated Projection Natural-Fit
Axle Plates:	Standard Amputee
Casters:	4", 5", 6", 8" composite 5", 6" soft-roll
Forks:	5" short fork 8" tall fork
Wheel Locks:	Pull-to-lock Push-to-lock Pull-to-lock with extensions Push-to-lock with extensions
Armrests:	Tubular swing-away (desk length) Height-adjustable drop-in (desk or full length)
Back Canes:	Straight, with or without push handles Angle-adjustable, with or without push handles 8 degree with push handles Straight reinforced
Backs:	Standard nylon sling with convenience pocket Tension-adjustable upholstery Synergy Shape Back
Front Riggings:	70 degree bi-directional swing-away leg rests Elevating leg rests with calf pads Comfort adjust footplates Nylon heel loops
Accessories:	Height-adjustable, removable anti-tip wheels Composite side guards Aluminum side guards Lap belt Nylon backpack Spoke protectors Transit securement package

1) Available via special order



(US) 866-800-2002 • (Canada) 888-570-1113 Visit us on the web: www.quantumrehab.com

The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.